

Software Name: Driver for CP-114EL Series (Windows)

Version No.	Build Date No./ Release Date	Supported OS	Apply to Models	New Features/ Enhancements	Bugs Fixed/Changes	Note
v1.21	Build 14012310	Windows_7, Windows_8, Windows_Server_2003, Windows_Server_2008, Windows_Server_2008_R2, Windows_Server_2012, Windows_XP, Windows_Vista, Windows_8.1, Windows_Server_2012_R2	PCIe/UPCI/PCI Serial Cards	New Features: 1. Supports Windows 8.1 2. Supports Windows Server 2012 R2 Enhancements: N/A	Bugs Fixed: 1. The serial mouse may not work when the system restores from sleep. 2. Fixed the issue that the driver cannot be installed due to the signature issue. Changes: N/A	N/A
v1.20	Build 12111310	Windows_2000, Windows_Server_2003, Windows_XP, Windows_Server_2008, Windows_Server_2008_R2, Windows_Vista, Windows_7, Windows_8, Windows_Server_2012	PCIe/UPCI/PCI Serial Cards	New Features: 1. Supports Windows 8 2. Supports Windows Server 2012 Enhancements: 1. Allow to use IOCTL_SERIAL_SET_FIFO_CONTR OL to clear RX FIFO (with value 0x02) or TX FIFO (with value 0x04.)	Bugs Fixed: 1. Flushing TX FIFO more times to prevent the serial port from sending out unknown data when the application closes the port. 2. Fixed linear baudrate so that it does not support below 200 bps. Changes: N/A	N/A